
WeightWatchers

POWERED BY

TANITA

Body Composition Monitor

Model: BC-730WW

en Instruction Manual

Read this Instruction Manual carefully and keep it for future reference.
Images in this manual are for illustrations purposes only.

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en Table of Contents

Introduction	1
Features and Functions.....	1
Preparations before Use	2
Getting Accurate Readings	2
Automatic Shut-down Function.....	2
Changing the Weight Mode	3
Setting (resetting) Date and Time	3
Setting and Storing Personal Data in Memory	4
Taking Weight and Body Fat Readings	5
Programming the Guest Mode	6
Taking Weight Measurement Only.....	6
Definition for Each Measurement	7
Troubleshooting	8

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Introduction

Welcome to your new Body Composition Monitor! Research shows that people who track their weight more frequently lose more weight.*

Remember, your weight is just one number that matters. That's why this scale measures body fat, body water and muscle mass—which, science says, are also key for overall health. It's important to also acknowledge (and celebrate!) your non-scale victories. Eating veggies? Clothes fitting looser? Have more energy? These wins all count.

The WeightWatchers programme is designed to help you reach and maintain your weight health goals. For more information visit ww.com today!

*Burke, L.A., et al. (2011). Self-monitoring in weight loss: a systematic review of the literature. Journal of the American Dietetic Association, 111(1), 92–102.

Note: Read this Instruction Manual carefully and keep it handy for future reference.

Safety Precautions

Do not use the Body Composition Monitor if you have an electronic medical implant such as a pacemaker, as it passes a low-level electrical signal through the body which may interfere with the operation of the implant.

Pregnant women should only use the weight function. Other functions are not intended for use when pregnant.

Do not use the product if there are any signs of damage (cracks, etc.) or missing parts on the exterior before or during use.

Important Notes for Users


The body fat monitor function is not intended for pregnant women, professional athletes or bodybuilders.

Recorded data may be lost if the unit is used incorrectly or is exposed to electrical power surges.

TANITA takes no responsibility for any kind of damage or loss caused by these units, or any kind of claim made by a third person.



Note: Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or overhydration due to such factors as alcohol consumption, menstruation, illness, intense exercise, etc.

Features and Functions


Accessories

AA-Size Batteries (4)

Preparations before Use

⚠ Inserting the Batteries

Open the battery cover on the back of the measuring platform. Insert the supplied AA-batteries as indicated.

Note: Be sure that the polarity of the batteries is set properly. If the batteries are incorrectly positioned the fluid may leak and damage floors. If you do not intend to use this unit for a long period of time, it is advisable to remove the batteries before storage.

Please note that the included batteries from factory may have decreased energy levels over time.

⚠ Positioning the Monitor

Place the measuring platform on a hard, flat surface where there is minimal vibration to ensure safe and accurate measurement.

Note: To avoid possible injury, do not step on the edge of the platform.

Handling Tips

This monitor is a precision instrument utilizing state-of-the-art technology.

To keep the unit in the best condition, follow these instructions carefully:

- Do not attempt to disassemble the measuring platform.
- Store the unit horizontally, and place it so that the buttons will not be pressed accidentally.
- Avoid excessive impact or vibration to the unit.
- Place the unit in an area free from direct sunlight, heating equipment, high humidity, or extreme temperature change.
- Never submerge in water. Use alcohol to clean the electrodes and glass cleaner (applied to a cloth first) to keep them shiny; avoid soaps.
- Do not step on the platform when wet.
- Do not drop any objects onto the platform.

Getting Accurate Readings

Heels centered on electrodes

Strategically placed electrodes ensure accurate measurements, accommodating foot sizes beyond the platform size.

To ensure accuracy, readings should be taken without clothing and under consistent conditions of hydration. If you do not undress, always remove your socks or stockings, and be sure the soles of your feet are clean before stepping on the measuring platform.

Be sure that your heels are correctly aligned with the electrodes on the measuring platform. Strategically placed electrodes ensure accurate measurements, accommodating foot sizes beyond the platform size.

It is best to take readings at the same time of day. Try to wait about three hours after rising, eating, or hard exercise before taking measurements.

While readings taken under other conditions may not have the same absolute values, they are accurate for determining the percentage of change as long as the readings are taken in a consistent manner. To monitor progress, compare weight and body fat percentage taken under the same conditions over a period of time.

Note: An accurate reading will not be possible if the soles of your feet are not clean, or if your knees are bent or you are in a sitting position.

Automatic Shut-down Function

The automatic shut-down function shuts off the power automatically in the following cases:

- If you interrupt the measurement process. The power will shut down automatically within 30 seconds depending upon the type of operation.
- During programming, if you do not touch any of the buttons within 60 seconds.*
- After you have completed the measuring process.

Changing the Weight Mode

You can change the weight unit by using the switch on the back of the scale.

Note: If weight mode is set to pounds or stone-pounds, the height programming mode will be automatically set to feet and inches. Similarly, if kilograms is selected, height will be automatically set to centimeters.

Setting (resetting) Date and Time

1. Set Date

Insert batteries into the scale and then press any button to start the Year setting. Press the <Up/Down> buttons to set year and then press SET button. The unit will beep once to confirm.


2. Set Hour

Press <Up/Down> buttons to set hour and then press SET button.
The unit will beep once to confirm.

3. Set Minute

Press <Up/Down> buttons to set minute and then press SET button.
The unit will beep twice to confirm the programming.
The power will then shut down automatically.

Note:

- After you replace the batteries or if you want to reset date and time, press and hold the Arrow (recall) button to begin the date and time programming process. Replacing the batteries will cause the loss of date and time settings.
- 60 seconds or more of inactivity will cause automatic shut off during this process. No data will be saved as a result.
- If you make a mistake or wish to turn the unit off before you complete the programming, press the OFF button to turn the unit off. No data will be saved as a result.

Setting and Storing Personal Data in Memory

The unit can be operated only if data has been programmed into one of the personal data memories.

1. Turn on the Power

Press the SET button to turn on the unit.

The unit will beep to confirm activation, the Personal data number will be displayed, and the display will flash.

Note: If you don't operate the unit for 60 seconds after the unit has been turned on, the unit will turn off automatically.

Note: If you make a mistake or want to turn the unit off before you have finished programming it, press the OFF button to force quit.

2. Select a Personal data Number


Press the <Up/Down> buttons to select a Personal data number (1, 2, 3, 4 and 5). Once you reach the Personal data number you wish to use, press the SET button. The unit will beep once to confirm.

3. Set your Birthday

Set the date of your birthday (Year / Month / Day) using <Up/Down> buttons. Press the SET button to confirm at each programming. The unit will beep once to confirm.

Note: The range of age is from 5–99.


4. Select Female or Male

Use the Up/Down buttons to scroll through Female (♀) or Male (♂), Female/Athlete (♀/🏃) and Male/Athlete (♂/🏃) settings, then press the SET button. The unit will beep once to confirm.

Who should use Athlete Mode?

Adults aged 18 years and over who either:


- Train or exercise for 12 hours or more a week and have been doing so for at least six months.
- Are body builders.
- Are professional athletes who want to monitor their progress at home.
- Have a lifetime history of fitness and who used to do more than 12 hours a week but do less now.

5. Specify the Height

The display defaults to 5' 70" (170cm) (range of height is from 3' 0.0"—7' 3.0" or 90cm–220cm).

Use the <Up/Down> buttons to specify Height and then press the SET button. The unit will beep once to confirm. The display shows the personal data one time to confirm, and then the display shows "0.0" and is ready for measurement.


6. Measuring your weight and body fat

Step on to the scale after "0.0" is shown on display within 30 seconds.

Do not step off until Body fat % is shown.

The readings will be shown for 40 seconds.

The unit will then shut off and the readings will be stored in the memory.


Note: If you step onto the scale before "0.0" appears the display will show "Err" and you will not obtain a reading.

Furthermore, if you do not step onto the measuring scale within about 30 seconds after "0.0" appears, the power is shut-off automatically. Please try the setting and storing personal data process again.

Note: If you change the birthday and/or gender, the unit will show "CLr".

If you select "YES" using <Up/Down> buttons, the existing personal data will be deleted.

If you select "no" using <Up/Down> buttons, no new data will be saved (existing data will be unchanged) then the unit will shut off.

Taking Weight and Body Fat Readings

After programming your personal data you are ready to take a reading.

1. Step onto the scale

Step onto the scale while the scale is off, the scale turns on automatically and starts measuring.


2. Get Your Readings

After the measurement is taken, readings will be automatically displayed in the order listed below. After the last result reading is displayed, the weight and Body Fat % reading will again appear, then the unit will automatically shut off.

The user may use the <Up> buttons to scroll thru the different results.

For advanced scrolling, the user may press and hold <Up> button.

Upon reading the desired readings, release the button to view the result, readings will appear for 40 seconds before the unit automatically shuts off.

If the wrong personal number appears on the display.

You can change the personal number by using the SET button.

Press the SET button until your personal number and profile appears.


After you select your correct personal number and wait for 4 seconds, the display will show your correct readings automatically.

The readings will be stored in your personal number memory.

If you select "dEL", no readings will be stored in memory.

Note: This product is not intended for people under the age of 18 years. For children (age 5–17), the unit will only display the Weight and the Body Fat % (it will not display the Healthy Range Indicator).

3. Recall Function

To obtain previous reading press Recall button while the readings are being displayed.

The following result of previous readings will be displayed by using the <Up> buttons.

- Weight
- Body Fat %
- Total Body Water %
- Muscle Mass
- Body Type
- Basal Metabolic Rate
- Metabolic Age
- Visceral Fat

Press Recall button again to go back to the current result display.

Note: This product is not intended for people under the age of 18 years. For children (age 5-17), the unit will only display the Weight and the Body Fat % (it will not display the Healthy Range Indicator).

Note: Do not step off until Body fat % is shown.

Programming the Guest Mode

Guest mode allows you to program the unit for a one-time use without reading a personal data number.

Press the Guest button to activate the unit. Then follow the steps: Set Age, Select Female or Male, and Specify Height.

Step onto the scale after "0.0" is shown on the display. Your weight will be shown first, continue standing on the scale until the body fat % reading is shown on the display. The reading will be shown for 40 seconds.

Note: Do not step off until Body fat % is shown.

Note: If you step onto the scale before "0.0" appears the display will show "Err" and you will not obtain readings.

Taking Weight Measurement Only

Press the Weight-only button. After 2 or 3 seconds, "0.0" will appear on the display and the unit will beep. When "0.0" appears on the display, step onto the scale. The body weight value will appear. The display will light for about 20 seconds, and then the power will automatically turn off.

Definition for Each Measurement

Weight	Weight is the total mass of your body in kilos / pounds.
Body Mass Index (BMI)	Body Mass Index (BMI) measures the ratio between your weight and height.
Body Fat	Body fat is measured in kilograms/pounds and is also shown as a percentage of your total body weight (body fat percentage). Body fat is essential for maintaining body temperature and for protecting your joints and internal organs. Energy is what your body needs, and this is supplied by food and drink. Energies are burned through physical activity and regular bodily processes. But if you consume excess energy, those excess energies will be stored in fat cells. If stored fat is not used as energy, it will create excess body fat. High levels of body fat can have a negative impact on your health in the long term. Strong evidence shows that having high levels of body fat increases the risk of high blood pressure, heart diseases, type 2 diabetes and certain types of cancer. Too little body fat can lead to osteoporosis later in life, as well as an irregular menstrual cycle and infertility in women. That's why it's important to monitor your body fat results with TANITA.
Muscle Mass	Muscle mass is the predicted muscle weight in your body, and includes skeletal muscles, cardiac and digestive muscles, and the water in these muscles. Muscles function as the engine to burn energy. When your muscle mass increases, the speed with which your body burns energy also increases. This, in turn, increases your BMR, and helps you to reduce excess body fat and lose weight healthily.
Body Type	This assesses muscle and body fat levels, and rates the result as one of nine types. As your activity level changes, the balance of body fat and muscle mass will gradually change, which in turn affects your overall physique.
Visceral Fat	Visceral fat is located deep in the core abdominal area, surrounding and protecting your vital organs. While having visceral fat is extremely important, having too much can increase the risk of heart disease, high blood pressure and type 2 diabetes. Even if your weight and body fat remain constant, as you get older, the distribution of fat changes, and it is more likely to shift to the abdominal area.
Basal Metabolic Rate (BMR)	BMR is the minimum amount of energy your body needs for vital functions such as your heart, lungs, and digestive system when rested (i.e. when not performing any physical activities). Your BMR is influenced by many factors such as your weight, muscle mass, age, and gender. BMR is a good starting point for your personal plan to move towards a healthier lifestyle.
Metabolic Age	The metabolic age is the result of the comparison between your BMR and your chronological age group. By exercising more you can improve your muscle growth, which will benefit your BMR.
Body Water	Body water is the percent of your weight that is made up of water. Body water can fluctuate throughout the day, week, and month for many different reasons, like exercise, how much fluid you drink, and/or menstrual cycles. While it is expected that body water will shift from day-to-day and week-to-week, the amount that it shifts should remain relatively consistent throughout your weight loss journey.

Troubleshooting

If the following problems occur ...then...

- **A wrong weight format appears followed by kg, lb, or st-lb**

Refer to "Changing the Weight Mode" on page 3.

- **"Lo" appears on the display, or all the data appears and immediately disappears.**

Batteries are low. When this message appears, be sure to replace the batteries immediately, since weak batteries will affect the accuracy of your measurements. Change all the batteries at the same time with new AA-size batteries.

Note: Your settings will not be erased from the memory when you remove the batteries.

- **"Err" appears while measuring.**

Please stand on the platform keeping movement to a minimum.

The unit cannot accurately measure your weight if it detects movement.

The Body Fat Percentage is more than 75%, readings cannot be obtained from the unit.

- **The Body Fat Percentage reading does not appear or "-----" appears after the weight is measured.**

Your personal data has not been programmed. Follow the steps on page 4.

Make sure socks or stockings are removed, and the soles of your feet are clean and properly aligned with the guides on the measuring platform.

- **"OL" appears while measuring.**

Reading cannot be obtained if the weight capacity is exceeded.

Specifications

BC-730WW

Weight Capacity

330lb (150kg) (23st 8lb)

Weight Increments

0.2lb (0.1kg) (0.2lb)

Applicable Age Range

Body Fat %

18–99 years old

Muscle Mass

18–99 years old

Body Type

18–99 years old

Visceral Fat

18–99 years old

Basal Metabolic Rate

18–99 years old

Metabolic Age

18–99 years old

Total Body Water %

18–99 years old

Personal Data

5 memories

Time accuracy

Within ± 1 minute per month (+18 – +28°C)

Power Supply

4 AA batteries (included)

Power Consumption

60 mA maximum

Disposal

This product is an electronic device. Please dispose of this product appropriately, not as general household waste. Be sure to follow the regulations in your area when disposing of this product.

Not allowed to mix batteries with consumer wastes!

As consumer you are legally bound to return used or discharged batteries. You can deposit your old batteries at the public collecting points in your town, or wherever the corresponding batteries are sold and specifically marked collecting boxes have been set up. In case of scrapping the apparatus, the batteries should be removed from it and deposited at the collecting points as well.

- Measurement categories other than body weight are values derived from factors such as electrical resistance in the body.
- For those exceeding the intended age group, the measurement results are for reference purposes only.
- The product design and specifications may be changed at any time without prior notice.

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*Burke, L.A., et al. (2011). Self-monitoring in weight loss: a systematic review of the literature.

Journal of the American Dietetic Association, 111 (1), 92–102.

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WW Logo and Weight Watchers are the trademarks of WW International, Inc. and are used under license by
© 2024 WW International, Inc. and TANITA Corporation All rights reserved.
EU Representative**TANITA Europe B.V.**
TANITA Corporation of America, Inc.
Arlington Heights, Illinois 60005 U.S.A.
TANITA Health Equipment H.K. Ltd.