WeightWatchers

POWERED BY

TANITA

Dual Frequency Body Composition Monitor

Model:RD-953WW

en

Instruction Manual

Read this Instruction Manual carefully and keep it for future reference.
Images in this manual are for illustrations purposes only.

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en Table of Contents

Introduction	3
Features and Functions	3
For Your Safety	4
Getting Accurate Readings	5
Preparations before Use	6
Connecting to Your Smartphone	7
Setting (Resetting) Date and Time	8
Setting and Storing Personal Data in the Memory	9
Taking Measurements	11
Measurement Results	13
Specifying the Personal Data Number for Measurement	14
Programming the Guest Mode	15
Measuring Weight Only	15
Definition for Each Measurement	16
Troubleshooting	17
Specifications	18

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Introduction

Welcome to your new Body Composition Monitor! Research shows that people who track their weight more frequently lose more weight.*

Remember, your weight is just one number that matters. That's why this scale measures body fat, body water and muscle mass—which, science says, are also key for overall health. It's important to also acknowledge (and celebrate!) your non-scale victories. Eating veggies? Clothes fitting looser? Have more energy? These wins all count.

The WeightWatchers programme is designed to help you reach and maintain your weight health goals. For more information visit ww.com today!

*Burke, L.A., et al. (2011). Self-monitoring in weight loss: a systematic review of the literature. Journal of the American Dietetic Association, 111(1), 92–102.

Note: Read this Instruction Manual carefully and keep it handy for future reference.

⚠ Safety Precautions

Do not use the Body Composition Monitor if you have an electronic medical implant such as a pacemaker, as it passes a low-level electrical signal through the body which may interfere with the operation of the implant.

Pregnant women should only use the weight function. Other functions are not intended for use when pregnant.

Important Notes for Users

The body fat monitor function is not intended for pregnant women, professional athletes or bodybuilders.



Recorded data may be lost if the unit is used incorrectly or is exposed to electrical power surges.

TANITA takes no responsibility for any kind of damage or loss caused by these units, or any kind of claim made by a third person.

Note: Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or overhydration due to such factors as alcohol consumption, menstruation, illness, intense exercise, etc.

Features and Functions


Front

Location Code

Back


Quick Start Guide

Accessories

AA-Size Batteries (4)

For Your Safety

This section explains precautionary measures to be taken to avoid injury to the users of this device and others, and to prevent damage to property. Please familiarise yourself with this information to ensure safe operation of this equipment.

	Warning	Failure to follow instructions highlighted with this mark could result in death or severe injury.
	Caution	Failure to follow instructions highlighted with this mark could result in injury or damage to property.
		This mark indicates actions that are prohibited.
		This mark indicates instructions that must always be followed.

Warning

 Prohibited	Do not use the product if there are any signs of damage (cracks, etc.) or missing parts on the exterior before or during use.
	This monitor must not be used by people with pacemakers or other medical implants. This monitor passes a weak electrical current through the body which could interfere with and cause the malfunction of electrical medical implants, resulting in serious harm to the user.
	Interpretation of measurements (e.g. evaluation of measurements and formulation of exercise programmes based on the measurements) must be performed by a professional. Implementing weight loss measures and exercise programmes based on self-analysis could be harmful to your health. Always follow the advice of a qualified professional.
	Do not step on the edge of the platform. Do not jump onto the monitor. Doing so may cause you to fall or slip, resulting in serious injury.
	Do not place this monitor on a slippery surface such as a wet floor. Doing so may cause you to fall or slip, resulting in serious injury. It may also cause water to get inside the monitor resulting in damage to internal components.
	Do not drop any objects onto the monitor. Avoid excessive impact or vibration to the monitor. Doing so may damage the glass, resulting in injury.

Caution

 Prohibited	Do not stand on the display screen or control panel.
	Do not insert your fingers into gaps or holes. Doing so may result in injury.
	Do not use rechargeable batteries. Do not use old batteries together with new batteries, or a mix of different types of batteries at the same time. Doing so may cause battery fluid to leak or the batteries to become excessively hot and rupture, resulting in damage to the monitor or injury.
	Never allow an infant or small child to carry the monitor. The child may drop the monitor, resulting in injury.
	This monitor is intended for home use only. This monitor is not intended for professional use including hospitals, medical or fitness facilities. It is not equipped with the components required for heavy usage. Using the monitor in a professional location will invalidate the warranty.

Place the monitor on a hard, flat surface where there is minimal vibration to ensure safe and accurate measurement.

Avoid excessive impact or vibration to the monitor.

Doing so may damage the monitor or cause it to malfunction, or may cause setting to be erased.

Do not use or store this monitor in the following locations.

- Locations that are exposed to direct sunlight, near to heaters, or exposed directly to air from air conditioning units, etc.
 - Locations where there are severe temperature changes or a large amount of vibration.
 - Locations with high levels of humidity or moisture.
- Doing so may damage the monitor or cause it to malfunction.

Do not step on the monitor when your body is wet (immediately after a shower or bath, etc.)

Be sure to thoroughly dry your body and feet before using the monitor. Failure to do so may cause you to slip or fall, resulting in injury. Also, the monitor cannot take accurate measurements if your body or feet are wet.

Keep away from water.

Do not use boiling water, benzene, thinner, etc. to clean the monitor.

If the monitor becomes dirty, soak a soft cloth in water or neutral household detergent, wring it out thoroughly and wipe the monitor clean, then use a dry cloth to wipe the monitor dry.

Do not attempt to disassemble the monitor. Do not modify this equipment in any way.

Do not wash the monitor in water.

Do not submerge this monitor in water.

Doing so may damage the monitor or cause it to malfunction.

Be sure to insert the batteries with the polarities (+, -) in the correct position.

If the batteries are incorrectly inserted, battery fluid may leak and damage the floor. (If you do not intend to use this monitor for a long period of time, it is advisable to remove the batteries before storage.)

Do not use the monitor near other products such as microwave ovens that emit electromagnetic waves.



Doing so may cause the monitor to malfunction, preventing accurate measurements.

Another person should assist persons with disabilities who may not be able to take a measurement alone.

Always...

Getting Accurate Readings

To ensure accuracy, readings should be taken without clothing and under consistent conditions of hydration. Always remove your socks or stockings, and be sure the soles of your feet are clean before stepping on the measuring platform. Be sure that the arch of each foot is aligned with the center of the measuring platform. Don't worry if your feet appear too large for the unit, accurate readings can still be obtained if your toes overhang the platform. It is best to take readings at the same time of day. Try to wait about two hours after getting up, eating, or strenuous exercise before taking measurements. While readings taken under other conditions may not have the same absolute values, they are accurate for determining the percentage change as long as the readings are taken in a consistent manner. To monitor your progress, compare your weight and body fat percentage readings taken under the same conditions over a period of time.


Note: An accurate reading will not be possible if the soles of your feet are not clean, if your knees are bent or if you are sitting down.

Preparations before Use

Inserting the Batteries

Note:


- Be sure that the batteries are inserted in the correct orientation. If the batteries are inserted incorrectly, fluid may leak from the batteries and damage the floor.
- If you do not intend to use this unit for a long period of time, remove the batteries before storage.
- The included batteries from factory may have decreased energy levels over time.

Set your location code

Different locations across the world have different gravitational pulls. These differences affect the accuracy of the weight measurement. You can obtain the most accurate weight readings by selecting the correct gravity setting on the scale according to your geographical location.

1. When the power is off, press and hold the Δ button.
2. Refer to the map enclosed, and press the $\wedge \vee$ buttons to select your location code.
3. Press the $\bigcirc$ button to enter and save the setting.

Changing the Weighing Unit

Note: Only for use with models that have the unit switching button.

The weighing unit can be changed using the switch on the back of the monitor.

Press the Δ button when the display is blank. When "0.00kg" (or "0.0lb", "0.0stlb") appears on the display, press the Unit Switching button underside of the monitor to switch the weighing unit.


Note: If weighing unit is set to pounds or stone-pounds, the height programming mode is automatically set to feet and inches. Similarly, if kilograms is selected, height is automatically set to centimeters.

Connecting to Your Smartphone

1. Download the MY TANITA Healthcare app in the Google Play store or the Apple App store. Also turn on **Bluetooth®** in the settings of your phone.

Note:

- If you are not using an app, enter settings manually in the Body Composition Monitor (page 9) before using.
- If you purchase a new smartphone, be sure to adjust the registration of your Body Composition Monitor accordingly.
- Approximate signal range: 5m in areas without obstructions

2. Launch the app, register with your email address and create a personal profile.


3. Click on CONNECT MY DEVICE and select your device.

4. Press and hold the O (communication) button on the device to start pairing. After pairing successfully, the app will display STEP ON. You can now start measuring.

Note:

The screen displays and instructions may differ depending on the app specifications.

5. Operate the Body Composition Monitor as indicated by the app. If the Body Composition Monitor is not powered on, press and hold the  button for 3 seconds or more, then enter your information as indicated by the app.


Note: Depending on signal strength, it may take a few seconds to connect to the Body Composition Monitor.

The  icon is lit while data is being transmitted.

To weigh yourself while holding a smartphone, first enter the amount of weight you want to subtract into "Weight of Clothing, etc." The input amount will be subtracted from your weight display.

If 0.20 kg is entered into "Weight of Clothing, etc.," the following will be displayed before a measurement is taken.

6. When the Body Composition Monitor is activated, "0.00kg" (0.0lb, 0.0stlb) is displayed.

7. After "0.00kg" (0.0lb, 0.0stlb) is displayed, carefully step barefoot onto the measuring platform.


Note: Please refrain from moving during measurement.

Registration of your body monitor is completed once measurement is completed.

Setting (Resetting) Date and Time (Not required if connection to the smartphone is complete)

1. Set Date

While the power is off, press the $\bigcirc$ button to start setting.
Press the $\wedge \vee$ buttons to set year and then press the $\bigcirc$ button.
Set the month and date in the same way.

Note:

- When the power is switched on for the first time, the date and time setting screen can also be displayed by pressing the $\wedge \vee$ buttons.
- When you have replaced the batteries, the date and time setting screen can also be displayed by pressing the $\wedge \vee$ buttons.

2. Set Time

Press the $\wedge \vee$ buttons to set hour and then press the $\bigcirc$ button.

Set the minutes in the same way.

The set date and time are displayed in order, and then the power automatically turns off.

Note:

- If the Body Composition Monitor is not operated for 60 seconds or more, the power will turn off automatically. Settings will not be saved in this case.
- If you make a mistake or wish to turn the monitor off before you complete the settings, press the $\triangle$ button to turn the monitor off. Settings will not be saved in this case.
- To change the date or time settings, press and hold the $\vee$ button for 3 seconds when the power is off, and then follow the procedure above to change the date or time settings.

Setting and Storing Personal Data in the Memory (Not required if connection to the smartphone is complete)

The Body Composition Monitor can only be operated if data has been programmed into one of the personal data memories. Only weight measurements can be taken if you use the Body Composition Monitor without registering personal data.

1. Select a Personal Data Number

Press the $\bigcirc$ button to turn on the monitor. Press the $\wedge \vee$ buttons to select a personal data number (1, 2, 3 or 4).

Press the $\bigcirc$ button to confirm the personal data number.

Note:

- The monitor turns off automatically if it is not operated for 60 seconds.
- If you make a mistake or want to turn the monitor off before you have finished programming it, press the $\triangle$ button to force quit. Partially registered data is saved on the monitor.
- The backlight color changes according to the personal data number.

- 1:Blue
- 2:Pink
- 3:White
- 4:Light blue

2. Set Your Birthday

Set the date of your birthday (Year / Month / Day) using the $\wedge \vee$ buttons.

Press the $\bigcirc$ button to enter each item.

Note:

Year settings can be configured from 1900 onwards.

3. Select Female or Male

Use the $\wedge \vee$ buttons to scroll through Male [♂], female [♀], male athlete [♂🏃] and female athlete [♀🏃] settings, then press the $\bigcirc$ button.

Who should use Athlete Mode?

Adults aged 18 years and over who either:

- Train or exercise for 12 hours or more a week and have been doing so for at least six months.
- Are body builders.
- Are professional athletes who want to monitor their progress at home.
- Have a lifetime history of fitness and who used to do more than 12 hours a week but do less now.

4. Set Your Height

The default height is 170cm (5' 7.0") (height range is 90cm - 220cm or 3'0.0" - 7'3.0").

Use the $\wedge \vee$ buttons to set your height and then press the $\bigcirc$ button.

The display shows "0.00kg" (or 0.0lb, 0.0stlb) and the monitor is ready for measurement.

Setting and Storing Personal Data in the Memory

5. Measuring Your Body Composition

Check "0.00kg" (0.0lb, 0.0stlb) is shown on the display, then step on the monitor.
Registered personal data and the measurement results are displayed in order, then the power switches off automatically.

Note:

If you step onto the monitor before 0.00kg is displayed, "Err,-OL" is displayed and the power turns off. Set the personal information again.

Also, if you do not step onto the monitor within 60 seconds of "0.00kg"(0.0lb, 0.0stlb) appearing, the power turns off automatically. In this case, you need to set and store personal data again.

Note:

The monitor needs to be connected to the smartphone in order to use the dedicated app.

Connect to the smartphone, and register the Body Composition Monitor as shown in the instructions on page 7.

If you change the birthday and/or gender, "CLr" is displayed on the monitor.

Select "YES" using the $\wedge \vee$ buttons and then press the $\bigcirc$ button, to delete the existing personal data enabling you to register new data.

If you select "no" using the $\wedge \vee$ buttons and then press $\bigcirc$ button, no new data is saved (existing data is left unchanged) and power turns off.

Taking Measurements

1. Go to MY MEASUREMENTS and click 'START MEASUREMENT'.
Follow instructions in the app.
When "0.00kg" (0.0lb, 0.0stlb) is displayed, step onto the monitor.


2. Measurement data is automatically sent to the smartphone.
If there is unsent data stored on the monitor, that data is also sent to the smartphone.

Note: When taking measurements while holding your smartphone, enter the weight that you want to subtract into the "Smartphone's weight" section in the dedicated app before taking measurements.
The weight that you entered to be subtracted from your weight is displayed.
The "Smartphone's weight" can be entered up to -2.0 kg (-4.4 lb, -4 lb).
The maximum weight (weighing capacity) is reduced by the amount entered for the "Smartphone's weight".

If you do not have your smartphone when taking measurements.

Simply step onto the monitor. The auto-recognition feature will automatically recall your personal data. Your monitor will store up to 10 sets of measurements.
Saved data is sent to the smartphone together with new measurements the next time the smartphone is connected to the monitor when taking measurements.

Note: Up to 10 measurements can be stored in the device per ID number.
Each new data entered after this amount causes the oldest data to be deleted.
🔔 will notify you if pairing is complete and there is unsent data.


Taking Measurements

You can start taking measurements when you have finished entering your personal data.

1. Step onto the Monitor

Step onto the monitor when the power is off. The monitor starts measuring automatically.

2. View Your Results

After taking measurements, your results are automatically displayed in the order listed on the next page.

After the last result data is displayed, the weight is displayed again and then the monitor automatically turns off.

Note:

- Do not step off until “0000” disappears.
- The backlight color changes according to the recognized personal data number. (Blue/Pink/White/Light blue)
- When the monitor is connected to the smartphone and measurement data is saved in the monitor,  lights in the personal data display and the measurement result display screen.
- If there are more than 6 sets of data stored, the number of saved sets of data is displayed after the personal data display lights.
- Up to 10 sets of measurement data can be stored in the monitor.
- If you have set a nickname in the dedicated app, the nickname is displayed instead of the personal data number when personal data is displayed.

If the wrong personal number is displayed

Press the  button to select your personal number while readings or personal data is being displayed.

Then press the  button to display your own personal data and measurement results.

Situations where mis-recognition may occur

- When measuring people who have similar weights and body shapes
 - When measurements have not been taken for a few days and there is a reasonable change in body weight
 - When measuring directly after having the monitor stood upright, carrying the monitor or storing it with items placed on top
- If the monitor continues to mis-recognize, specify your personal data number manually before measurement. (page 14)

Measurement Results

Weight

Body Mass Index(BMI)

Body Fat %

Muscle Mass

**Body Type
(Physique Rating)**

Visceral Fat

Basal Metabolic Rate

Metabolic Age

Total Body Water %

Note:

- This product is not intended for people under the age of 18 years.
- For children (age 5-17), the monitor only displays weight, BMI, and body fat percentage.
- For children (age 4 or less) and adults (age over 100), the monitor only displays the weight.
- "Err" is displayed if the body fat percentage exceeds the measurable range.
- You can use the  button to scroll through the results.
- Press and hold the  button for fast scrolling.
- After reading the results that you want to see, release the button to see the rest of the results.
- The monitor automatically turns off after 40 seconds.
- Press the  button to cancel measurement or turn the power off.
- Do not move the monitor for approximately 3 seconds after taking measurements and stepping off.
To prevent accidents such as falling off the monitor, always step off the monitor before operating the buttons.

Specifying the Personal Data Number for Measurement

Follow these instructions if the monitor continues to mis-recognize your number.

1. When the power is off, press the $\wedge$ button to start measuring.
Select your personal number using the $\wedge$ / $\vee$ buttons and then press the $\bigcirc$ button.

2. Step onto the monitor when "0.00kg" (0.0lb, 0.0stlb) appears in the display.
The monitor turns off automatically after displaying the measurement results.

Note:

- When the monitor is connected to the smartphone and measurement data is saved in the monitor,  lights in the personal data display screen and the measurement result display screen.
- If there are more than 6 sets of data stored, the number of saved sets of data is displayed after the personal data display lights.
- Up to 10 sets of measurement data can be stored in the monitor.
- If you step onto the monitor before 0.00kg (0.0lb, 0.0stlb) is displayed, "Err,-OL" is displayed and the power turns off during measurement.

Programming the Guest Mode

Guest mode allows you to program the monitor for one-time use without using a personal data number.

Press the $\wedge$ button to select "GUEST", then press the $\bigcirc$ button.

Next, follow the steps below:

Set the age, gender (male/ female/ male athlete/ female athlete) and height.

Step onto the monitor after "0.00kg" (0.0lb, 0.0stlb) is displayed.

The readings are displayed for 40 seconds.

Note:

The age range is from 5 to 99.

Who should use Athlete Mode?

Adults aged 18 years and over who either:


- Train or exercise for 12 hours or more a week and have been doing so for at least six months.
- Are body builders.
- Are professional athletes who want to monitor their progress at home.
- Have a lifetime history of fitness and who used to do more than 12 hours a week but do less now.
- Do not step off until "0000" disappears.

Measuring Weight Only

Press the $\triangle$ button. After 2 or 3 seconds, "0.00kg" (0.0lb, 0.0stlb) is displayed.

When "0.00kg" (0.0lb, 0.0stlb) is displayed, step onto the monitor.

The body weight measurement is displayed, and the monitor power turns off automatically after approximately 20 seconds.

Definition for Each Measurement

Weight	Weight is the total mass of your body in kilos / pounds.
Body Mass Index (BMI)	Body Mass Index (BMI) measures the ratio between your weight and height.
Body Fat	Body fat is measured in kilograms/pounds and is also shown as a percentage of your total body weight (body fat percentage). Body fat is essential for maintaining body temperature and for protecting your joints and internal organs. Energy is what your body needs, and this is supplied by food and drink. Energies are burned through physical activity and regular bodily processes. But if you consume excess energy, those excess energies will be stored in fat cells. If stored fat is not used as energy, it will create excess body fat. High levels of body fat can have a negative impact on your health in the long term. Strong evidence shows that having high levels of body fat increases the risk of high blood pressure, heart diseases, type 2 diabetes and certain types of cancer. Too little body fat can lead to osteoporosis later in life, as well as an irregular menstrual cycle and infertility in women. That's why it's important to monitor your body fat results with TANITA.
Muscle Mass	Muscle mass is the predicted muscle weight in your body, and includes skeletal muscles, cardiac and digestive muscles, and the water in these muscles. Muscles function as the engine to burn energy. When your muscle mass increases, the speed with which your body burns energy also increases. This, in turn, increases your BMR, and helps you to reduce excess body fat and lose weight healthily.
Body Type (Physique Rating)	This assesses muscle and body fat levels, and rates the result as one of nine types. As your activity level changes, the balance of body fat and muscle mass will gradually change, which in turn affects your overall physique.
Visceral Fat	Visceral fat is located deep in the core abdominal area, surrounding and protecting your vital organs. While having visceral fat is extremely important, having too much can increase the risk of heart disease, high blood pressure and type 2 diabetes. Even if your weight and body fat remain constant, as you get older, the distribution of fat changes, and it is more likely to shift to the abdominal area.
Basal Metabolic Rate (BMR)	BMR is the minimum amount of energy your body needs for vital functions such as your heart, lungs, and digestive system when rested (i.e. when not performing any physical activities). Your BMR is influenced by many factors such as your weight, muscle mass, age, and gender. BMR is a good starting point for your personal plan to move towards a healthier lifestyle.
Metabolic Age	The metabolic age is the result of the comparison between your BMR and your chronological age group. By exercising more you can improve your muscle growth, which will benefit your BMR.
Body Water	Body water is the percent of your weight that is made up of water. Body water can fluctuate throughout the day, week, and month for many different reasons, like exercise, how much fluid you drink, and/or menstrual cycles. While it is expected that body water will shift from day-to-day and week-to-week, the amount that it shifts should remain relatively consistent throughout your weight loss journey.

Troubleshooting

If the following problems occur, follow the instructions below.

"Lo" appears on the display.

Check the orientation of the batteries, and insert them correctly.

The batteries are running low. Replace the batteries promptly.

Refer to page 6.

"OL" appears on the display during use

The maximum weight (weighing capacity) 200kg is exceeded. (Measurements cannot be taken)

"Err" appears on the display after measurement

Are the personal data settings correct?

Did you select the wrong personal data number?

Measurement items have exceeded the measurement range. (Measurements cannot be taken)

Cannot measure body composition even though personal data is set

The personal data settings are not complete. You need to measure body composition once immediately after entering the personal data.

Cannot connect my smartphone to the Body Composition Monitor

"Err pair" or "Err DATA" appears on the display

Are batteries inserted correctly in the Body Composition Monitor?

Are the batteries running low?

Are the smartphone and body composition meter within communication range?

The communication range for The Body Composition Monitor is 5m in an unobstructed location.

Is the "Bluetooth" setting in the smartphone "Settings" turned on?

"Err pair" appears repeatedly.

Delete the pairing information from "Settings" > "Bluetooth" on the smartphone, and pair the smartphone with the Body Composition Monitor again.

"Err UUID" appears repeatedly.

Delete the pairing information from "Settings" > "Bluetooth" on the smartphone, and pair the smartphone with the Body Composition Monitor again.

"Err c9" appears on the display

Remove and replace the batteries from the body composition meter, and set the date and time again.

The date and time settings have been changed

When using the Body Composition Monitor linked with an smartphone, the date and time settings on the Body Composition Monitor are overwritten with those from the smartphone.

The "Input Data" button in the app does not work

Is the "Bluetooth" setting in the smartphone "Settings" turned on?

Linking between the app and the Body Composition Monitor is only enabled when the Body Composition Monitor is turned off.

Check that the power of the Body Composition Monitor is turned off and try again.

Specifications

RD-953WW

Weight Capacity

200 kg (440 lb)(31st 6 lb)

Weight Increments0-100 kg/0.05 kg 100-200 kg/0.1 kg
0-200 lb/0.1 lb 200-440 lb/0.2 lb 0-31st 6 lb/1 lb**Body Fat %**

18-99 years old

Body Mass Index

18-99 years old

Muscle Mass

18-99 years old

Body Type (Physique Rating)

18-99 years old

Visceral Fat

18-99 years old

Basal Metabolic Rate

18-99 years old

Metabolic Age

18-99 years old

Body Water %

18-99 years old

Personal Data

4 memories

Time accuracyWithin ± 1 minute per month (+18 – +28°C)**Power Supply**

4×AA Alkaline batteries (included)

Power Consumption

250mA maximum

Communication Method

Bluetooth Low Energy technology

Frequency Band

2402 – 2480 MHz

Maximum RF Output Power (adaptive frequency)

2 dBm (2402-2480 MHz)

Disposal

This product is an electronic device.

Please dispose of this product appropriately, not as general household waste. Be sure to follow the regulations in your area when disposing of this product.

⚠ Not allowed to mix batteries with consumer wastes!

As consumer you are legally bound to return used or discharged batteries. You can deposit your old batteries at the public collecting points in your town, or wherever the corresponding batteries are sold and specifically marked collecting boxes have been set up. In case of scrapping the apparatus, the batteries should be removed from it and deposited at the collecting points as well.

Hereby, TANITA declares that the radio equipment type RD-953WW is in compliance with Directives: 2014/53/EU and 2011/65/EU.

The full text of the EU declaration of conformity is available at the following internet address:

<https://tanita.eu/declaration-of-conformity>

- Measurement categories other than body weight are values derived from factors such as electrical resistance in the body. BMI is a calculated value.
- For those exceeding the intended age group, the measurement results are for reference purposes only.
- The product design and specifications may be changed at any time without prior notice.
- Apple, the Apple logo, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.
- Android and Google Play are trademarks of Google LLC.
- The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by TANITA Corporation is under license.

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*Burke, L.A., et al. (2011). Self-monitoring in weight loss: a systematic review of the literature.

Journal of the American Dietetic Association, 111 (1), 92–102.

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WW Logo and Weight Watchers are the trademarks of WW International, Inc. and are used under license by
© 2024 WW International, Inc. and TANITA Corporation All rights reserved.
TANITA Europe B.V.**
TANITA (Shanghai) Trading Co., Ltd.
TANITA Corporation of America, Inc.
U.S.A.